EN - English
SV - Svenska
DA - Danske
NO - Norsk
FI - Suomi

GENERAL SAFETY INSTRUCTIONS

READ ALL INSTRUCTIONS BEFORE USE

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE

When using electrical products, basic safety precautions should always be followed, including the following:

EN

1. The appliance is only intended for domestic/private use, not for commercial use.
2. Never use the appliance if damaged or if any of its accessories are damaged.
3. Never use the appliance during a thunderstorm.
4. In the event of defect or malfunction, switch off immediately and disconnect from the power supply and the mains outlet.
5. Never pull the power cable when disconnecting from the power supply or mains outlet.
6. Always pull the plug. Never hold or carry the appliance by the power cable.
7. Keep cables away from hot surfaces.
8. Always unplug the appliance from the electrical outlet after using and before maintenance. Turn all controls to the "OFF" position and remove the plug from the outlet.
9. Never leave the appliance unattended when plugged in. Unplug from the mains outlet when not in use and before changing or taking off parts or attachments.
10. The appliance can be used by persons from 15 years of age and above and persons with impaired physical, sensory or mental capabilities or lack of experience and knowledge if supervised or instructed in safe use of the appliance and they understand the hazards involved.
11. Children must not play with the appliance.
12. Do not touch an appliance that has fallen in water or other liquids. Switch off at the mains and unplug immediately. Keep dry – do not operate in wet or moist conditions.
13. Do not insert fingers, pins, needles, or sharp objects.
14. Do not open the appliance, even using the zip fastener (fitted for reasons relating to manufacture).
15. Make sure all moving parts can move freely at all times.
16. Do not sit or stand on the moving parts with all your weight as it may damage the appliance.
17. A damaged cable must be replaced by authorized service personnel only.
18. Never use under covers such as pillows or blankets.
19. Never use close to any flammable materials.
20. Keep the appliance away from high temperatures.

SPECIAL SAFETY PRECAUTIONS

1. Unsuitable users: The following users must not use the appliance or consult their doctor before use: Anyone under medical treatment or feeling unfit, suffering from hypertension, heart diseases, brain diseases or thrombosis, malignant tumors, pregnant, suffering from osteoporosis, with implanted pacemaker or other medical device.
2. Do not operate without carefully reading this operating manual.
3. Do not allow children or pets to play around the massager to avoid injury during operation of the massager.
4. Always wait 20 minutes before reusing the appliance.
5. Never use when feeling unwell, if parts are damaged, if smells or smoke occur during use, during a thunderstorm.
6. Do not use while driving.
7. Do not use if feeling drowsy and do not fall asleep during use.
8. Do not connect the massager to any other power supply or other electrical equipment at the same time.
9. Keep the massager away from water or liquids.
10. Do not place the massage heads directly on your hand, elbow, knee, or chest.
11. CAUTION: Do not use this product in dusty, hot places

Parts

1. Power Button
2. Speed control Button
3. 2 different massage heads

—3

EN

Operating instructions

1. Check the massager is secure and safe to use.
2. Insert the power plug into a proper wall outlet.
3. Press ON/OFF button to start operating with default low speed.
4. Press M button to switch into the high speed; press M button again to switch into low speed.
5. Press ON/OFF button to switch off the massager.
6. Switch off power and remove power plug from wall outlet.
7. The appliance may be operated for 15 minutes, before a pause of 20 minutes is recommended.

Care & maintenance

- Disconnect the power supply from the mains socket and the appliance before cleaning.
- Clean using a moist cloth or sponge only with soft detergent for delicate fabrics. Never use cleaning agents containing solvents.
- Make sure no fluid enters the appliance or accessories when cleaning.
- Never use the appliance until it is completely dry.

Troubleshooting

Problem	Cause	Solution
The item does not work	1. The product has just finished 20 minutes of massage, and automatically turn off. 2. The product has been overloaded. Too much weight on the massage heads.	1. Press on/off button to start again. 2. Turn off remove the weight, unplug the power cord, wait until cool down, turn on again.
The item is too hot	1. The heating function is on.	1. Turn off the heating button.

Warranty

Champion Nordic guarantees that this product is free from manufacturing defects in terms of materials and workmanship during a period of 2 years in Sweden, Norway, Denmark and Finland. This warranty period begins on the day the product is purchased or delivered. Champion Nordic has no obligation to repair or replace products that are not accompanied by a valid proof of purchase. This warranty applies only to products purchased and used for home use and does not cover damage resulting from abuse, failure to comply with the instructions from Champion Nordic, cases in which the product has been modified or subjected to unauthorized repair, improper packaging on the part of the owner, normal wear and tear, or improper handling on the part of a transport company.

Correct Disposal of this product

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources.

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